

HAYES

THE LOCAL HANDPIECE REPAIR CENTER

OneStep™ Handpiece Lubricant



NEW!
Metered-Dose*

NO MESS

Our superior patent pending metered-dose system means you'll use the right amount of lubricant every time

SAVE MONEY

Over 3,000 spritzes in each can! Just 1 PENNY per use!

PROTECTS HANDPIECES

Received top rating on industry-standard 4-ball wear test, which means longer turbine life

Also available:

- dropper bottles
- standard aerosol
- 500 ml size for Adec Assistina

2-PACK
JUST \$78

#HA-30.04


Clean AND lubricate in 1 easy step!

*patent pending


OneStep™ Handpiece Lubricant

Choose the Right Adapter

Our Metered-Dose Aerosol is
perfect for all handpieces!



Star—#HA-30.06



NSK—#HA-30.11



Lares—#HA-30.07



Kavo—#HA-30.13



W&H/Adec less
bulb—#HA-30.08



Midwest XGT/Stylus/
Bien-Air—#HA-30.14



W&H/Adec w/bulb/
Schein—#HA-30.09



5-Hole—#HA-30.16



E-Type attach—#HA-30.10



3-Hole—#HA-30.17

Simply thread the appropriate adapter on
the white cap and lubricate as directed



www.hayeshandpiece.com